



**Application Notes:** For application of Cas9 for site-specific genome editing in eukaryotic cells and organisms, we have computationally identified suitable target sites for the *S. pyogenes* Cas9. These sites are viewable as [UCSC Genome Browser tracks](#) for the human, mouse, rat, zebrafish, *C. elegans*, and *D. melanogaster* genomes. Sites are selected such that the seed sequence for each SpCas9 target site, 5'-NNNNNNNNNNNN-NGG-3', is specific to the relevant genome. A protocol for oligo cloning is included here and the plasmid sequences and additional information can be found at Zhang Lab's website ([www.genome-engineering.org](http://www.genome-engineering.org)).

**Citation Information:** please reference the following publication for the use of these material.

*Multiplex Genome Engineering using CRISPR/Cas Systems*


Cong L, Ran FA, Cox D, Lin S, Barretto R, Habib N, Hsu PD, Wu X, Jiang W, Marraffini LA, Zhang F. Science. 2013 Jan 3. DOI: 10.1126/science.1231143

## Target Sequence Cloning Protocol

(standard de-salted oligos are sufficient)


### pX260 (or pX334) – hSpCas9 (or hSpas9n nickase) + CRISPR array + tracrRNA:

In order to clone the target sequence into the pX260 backbone, synthesize two oligos of the form:



### pX330 (or pX335) – hSpCas9 (or hSpCas9n nickase) + chimeric guideRNA:

In order to clone the target sequence into the pX260 backbone, synthesize two oligos of the form:



\* \* \* \* \*

### Oligo annealing and cloning into backbone vectors (original version, new version on next page):

1. Digest 1ug of pX260 or pX330 with *BbsI* for 30 min at 37C:

1 ug	pX260 or pX330
1 ul	FastDigest <i>BbsI</i> (Thermo)
1 ul	FastAP (Thermo)
2 ul	10X FastDigest Buffer
X ul	ddH <sub>2</sub> O
<hr/>	
20 ul	total

2. Gel purify digested pX260 or pX330 using QIAquick Gel Extraction Kit and elute in EB.

3. Phosphorylate and anneal each pair of oligos:

1 ul	oligo 1 (100μM)
1 ul	oligo 2 (100μM)
1 ul	10X T4 Ligation Buffer (NEB)
6.5 ul	ddH <sub>2</sub> O
0.5 ul	T4 PNK (NEB)
<hr/>	
10 ul	total

Anneal in a thermocycler using the following parameters:

37°C	30 min
95°C	5 min and then ramp down to 25°C at 5°C/min

4. Set up ligation reaction and incubate at room temperature for 10 min:

X ul	<i>BbsI</i> digested pX260 or pX330 from <b>step 2</b> (50ng)
1 ul	phosphorylated and annealed oligo duplex from <b>step 3</b> (1:250 dilution)
5 ul	2X Quickligation Buffer (NEB)
X ul	ddH <sub>2</sub> O
<hr/>	
10 ul	subtotal
1 ul	Quick Ligase (NEB)
<hr/>	
11 ul	total

5. (optional but highly recommended) Treat ligation reaction with PlasmidSafe exonuclease to prevent unwanted recombination products:

11 ul	ligation reaction from <b>step 4</b>
1.5 ul	10X PlasmidSafe Buffer
1.5 ul	10mM ATP
1 ul	PlasmidSafe exonuclease
<hr/>	
15 ul	total

Incubate reaction at 37C for 30 min.

6. Transformation

## Oligo annealing and cloning into backbone vectors (new version with single-step digestion-ligation)

The following protocol is a new version that is simpler to set up and increase cloning efficiency (higher percentage of correct colonies).

**[IMPORTANT NOTE] Due to the simultaneous digestion-ligation step, the guide oligos CANNOT contain any *BbsI* enzyme site (i.e. the nucleotide sequence 'GAAGAC' or 'GTCTTC'). Please double-check your guides!**

1. Phosphorylate and anneal each pair of oligos:

1 ul	oligo 1 (100 $\mu$ M)
1 ul	oligo 2 (100 $\mu$ M)
1 ul	10X T4 Ligation Buffer (NEB)
6.5 ul	ddH <sub>2</sub> O
0.5 ul	T4 PNK (NEB)
10 ul	total

Anneal in a thermocycler using the following parameters:

37°C	30 min
95°C	5 min and then ramp down to 25°C at 5°C/min

Dilute the annealed oligo 1:250 (250-fold).

2. Set up digestion-ligation reaction:

X ul	pX330 or other backbone vector (100ng)
2 ul	phosphorylated and annealed oligo duplex from <b>step 1</b> (1:250 dilution)
2 ul	10X Tango buffer (or FastDigest Buffer)
1 ul	DTT (10mM to a final concentration of 1mM)
1 ul	ATP (10mM to a final concentration of 1mM)
1 ul	FastDigest <i>BbsI</i> (Thermo Fisher Fermentas)
0.5 ul	T7 DNA ligase
Y ul	ddH <sub>2</sub> O
20 ul	total

Incubate the ligation reaction in a thermocycler:

37°C	5 min
23°C	5 min
	Cycle the previous two steps for 6 cycles (total run time 1h)
4°C	hold until ready to proceed

3. (optional but highly recommended) Treat ligation reaction with PlasmidSafe exonuclease to prevent unwanted recombination products:

11 ul	ligation reaction from <b>step 4</b>
1.5 ul	10X PlasmidSafe Buffer
1.5 ul	10mM ATP
1 ul	PlasmidSafe exonuclease
15 ul	total

Incubate reaction at 37C for 30 min.

4. Transformation with 1-2 ul of the final product into competent cells

5. Pick colony and sequence verify the clones.